NEXUS™  | *Aortic Arch Stent
Graft System*

INSTRUCTIONS FOR USE

TABLE OF CONTENTS

| | | |
|----|---|-----------|
| 1 | Abbreviations | 2 |
| 2 | Device description | 3 |
| 3 | Indication for use | 9 |
| 4 | Contraindications | 9 |
| 5 | Warnings and precautions | 9 |
| 6 | Anticipated adverse events | 12 |
| 7 | Treatment and follow-up | 13 |
| 8 | How supplied | 14 |
| 9 | Clinical use information | 14 |
| 10 | Directions for use | 15 |
| 11 | Bail-out techniques | 23 |
| 12 | Follow-up imaging recommendation | 24 |
| 13 | Explanation of symbols | 24 |
| 14 | Appendix A – Consumables recommendation | 25 |
| 15 | Disclaimer of warranty | 26 |

1 | ABBREVIATIONS

| | |
|----------------------------|--|
| AAA | Abdominal Aortic Aneurysm |
| ARDS | Acute Respiratory Distress Syndrome |
| ASC | Ascending |
| BCA | Brachiocephalic Artery |
| BCT | Brachiocephalic Trunk |
| Caudal | Toward the legs |
| Cranial | Toward the head |
| CVA | Cerebral Vascular Accident |
| DE | Descending Extension |
| DESC | Descending |
| Distal (delivery system) | Distant, relatively to the user |
| Distal (implant) | Distant, relatively to the patient's heart |
| DS | Delivery System |
| EtO | Ethylene Oxide |
| Fr | French (units) |
| IFU | Instruction For Use |
| IMH | Intramural Hematoma |
| LAO | Left Anterior Oblique |
| LCCA | Left Common Carotid Artery |
| LSA | Left Subclavian Artery |
| MI | Myocardial Infarction |
| MOF | Multiple organ failure syndrome |
| MRI | Magnetic Resonance Imaging |
| NiTi | Nitinol |
| Proximal (delivery system) | Close, relatively to the user |
| Proximal (implant) | Close, relatively to the patient's heart |
| RAO | Right Anterior Oblique |
| RCCA | Right Common Carotid Artery |
| RHP | Rapid Heart Pacing |
| RSA | Right Subclavian Artery |
| SAR | Specific Absorption Rate |
| SG | Stent Graft |
| SPS | Self-Projecting Sleeve |
| STJ | Sino-Tubular Junction |
| TA | Tantalum |
| TIA | Transient Ischemic Attack |
| WBA | Whole Body Average |

2 | DEVICE DESCRIPTION

The Nexus™ Aortic Arch Stent Graft System, named hereinafter, Nexus™, is indicated for the endovascular treatment of thoracic aortic diseases involving the aortic arch.

The Nexus™ system excludes the diseased area and provides a permanent alternative conduit for arterial blood flow.

2.1 Nexus™ Stent Graft

The Nexus™ aortic arch stent graft system is comprised of two components of implantable stent grafts both in disposable 20 Fr. delivery system.

The implantable stent graft system is composed of two elementary and optional multiple extension stent graft modules, each introduced separately into the patient's vascular system:

- The Main Module, whose cranial narrow end is intended to be deployed into the Brachiocephalic artery and whose distal end is to be deployed into the Descending Thoracic Aorta.
- The Ascending Module, intended to be deployed in the Ascending Aorta.

The pre-loaded stent graft components are sequentially advanced to the diseased location over a guide wire using fluoroscopic guidance. Main module first and ascending to follow. Upon deployment, the stent graft components are collectively intended to conform to the shape and size of the proximal and distal seal zones surrounding the diseased location (**Figure 1**).


Figure 1: Configuration of the Nexus™ Main & ASC Stent Grafts within the anatomy

- Optional - Descending Extension can be used in case the aortic lesion elongates further distally and out of the covered length offered by the Main Module.

Each stent graft module is made of Nickel-Titanium (Nitinol) alloy stents, sewn to a polyester fabric, using surgical ePTFE coated polyester suture material. According to the clinical outcome assessment, both B-markers and Dot markers are made of tantalum (>99.9% purity) with no trace of Platinum. Radiopaque markers aid in visualization of the stent graft under fluoroscopy to facilitate accurate placement of the device components relative to each other and in relation to the patient’s vascular system. Fabric is inside the stent skeleton.

Table 1 – Stent Graft Materials

| COMPONENT | MATERIAL |
|------------------------------------|--|
| Stents | Nickel-Titanium (Nitinol) Alloy |
| Button Radiopaque and Ring Markers | Tantalum |
| “B” Markers | Tantalum |
| Graft Material | Polyester (PET= Polyethyleneterephthalate) |
| Suture | Polyethylene Polytetrafluoroethylene (PTFE) |

The stent graft is designed to be placed in the native vessel such that the unconstrained stent graft diameter is larger than the internal diameter of the native vessel. This “oversizing” helps excluding the lesion from the aortic blood flow and ensures that the stent graft is held in place.

The range of diameters and lengths is described in Section 2.2 Device Size, Table 2 – Graft Sizing Chart.

2.1.1 Main Module (Brachiocephalic Trunk to Descending Aorta)

The Main Module is the base stent graft, implanted from the BrachiocephalicTrunk (cranially) to the descending Aorta (caudally), and is aligned with the Aortic arch region (Figure 2). The shape of the module within the body is a thunderbolt shape (⚡) which contributes to its stability and fixation. This module contains a sleeve named the “Self-projecting sleeve” (SPS) with a single fenestration which opens toward the ascending Aorta.


Figure 2: Main module (Brachiocephalic to Descending Aorta)

The positioning of the module should be such to allow blood flow to the **Brachiocephalic** artery and to the descending aorta. Blood supply to the Left carotid artery should be done through bypass. Blood supply to the Left subclavian artery should be done through bypass or it can be sacrificed, according to physicians discretion. The Main Module implantation is done over a Bracheo-femoral wire (known as “Through & Through” technique) in which a stiff guide wire is placed into the vasculature from the Axillary/Brachial artery to the iliac-femoral of the patient (**Figure 3**).


Figure 3: Through & Through Technique

2.1.2 Ascending Aorta Module

The Ascending Module is the most proximal module of the stent graft system, which is placed in the Ascending Aorta, distally to the Coronary Arteries and the Sino-Tubular Junction (STJ). There are 2 configurations for the ascending module, a straight configuration and an oriented configuration.

Note: The straight ascending module configuration is preferable in anatomies where the proximal landing area is straight (for example if the ascending was previously replaced or wrapped), whereas the oriented ascending module better accommodates the curved proximal landing area. Choosing between straight and oriented ascending module is according to surgeon’s discretion.

The proximal strut of the Ascending Module is slightly inwardly bent to prevent damage to the Ascending Aorta (**Figure 4**). The distal strut of the module is partially exposed and contains non traumatic, triangular shaped fixation elements. These fixation elements are intermittently proximally and distally oriented, so as to provide bidirectional locking of the Ascending Module inside the SPS of the Main Module. Once deployed it engages the Main Module, thereby achieving a secure connection between the two modules.


Figure 4: Ascending Module

2.1.3 Descending Extension

The Descending Extension (DE) stent graft module is used in cases where the aortic lesion elongates further distally and out of the covered length offered by the Main Module stent graft. The DE stent graft module can be used in continuation to the Main Module to ensure exclusion of the lesion from the blood flow. Multiple DE modules may be utilized.

The DE is comprised of Polyester fabric on which equally spaced struts are sewn together with Tantalum radiopaque markers by two types of medical grade threads – the last two struts at each end are sewn from inside the graft, two opposite dot markers indicate both ends of the prosthesis and a “B” shaped marker is positioned at the minimally required overlap length. In all DE modules at both ends of the prosthesis the graft follows the end strut’s exact circumference to create a crown-like shape.

The DE module is available in two configurations: Cylindrical and Tapered (Figure 5). Please refer to section 2.2 for device sizes.


Figure 5: Both configurations of the Descending Extension Module

2.2 Device size

The Nexus™ stent graft systems are available in the sizes described in the Table below.



For questions about device sizing, refer to contact information in the back of this Instructions for Use, info@endospan.com.

Stent Graft Sizing Chart

| Sizing Chart Main Module: Brachiocephalic (BCA) to Descending Aorta (DESC) | | Sizing Chart – Ascending (ASC) Module: | | |
|---|-------------|--|----------------------|----------------------|
| Parameter | Size (mm) | Parameter | Straight - Size (cm) | Oriented - Size (cm) |
| \varnothing_{BCA} | 14/17/20 | \varnothing_{ASC} | 36/40/43 | 36/40/43 |
| L_{BCA} | 20/30/40 | L nominal | 50*/60/75 | 40*/55/70 |
| \varnothing_{DESC} | 32/36/40/44 | | | |
| L_{DESC} | 160/180 | | | |
| | | | | |

Note: All module sizes are compatible with each other.

* For the shortest ascending module, ascending Aorta landing zone diameter should be: $29 \leq \varnothing \leq 36$ mm.

| Sizing Chart – Descending Extension (DE): Cylindrical Module | | Sizing Chart – Descending Extension (DE): Tapered Module | |
|---|-------------|---|-------------|
| Ø Diameter (mm) | Length (mm) | Ø Diameter (PROX X DIST) | Length (mm) |
| 31 | 125/157/189 | 36 X 31 | 125/157/189 |
| 36 | | 40 X 36 | 125/157/189 |
| 40 | | 43 X 40 | 130/162/194 |
| 43 | 132/166/200 |  | |
|  | | | |

2.3 Nexus™ Delivery System

The Nexus™ Delivery System, which delivers all stent graft configurations, consists of a single-use, disposable hydrophilically coated catheter with an integrated handle to provide controlled deployment. It is available in 20 Fr. (outer diameter) for all sizes, which enables a percutaneous approach. The catheter assembly is flexible and compatible with a 0.035 inch (0.89 mm) guide wire. **(Figure 6).**


Figure 6 –Nexus™ Delivery System Handle

The Nexus™ delivery system (**Figure 7.1**) includes a set of 2 single-use disposable catheters, one for each module:

1. Main module delivery System
2. Ascending module delivery System
3. Optional- if extended coverage is required for the descending aorta, a third delivery system containing the Descending Extension Module should be made available.

The configuration of the main module delivery system is pre-shaped (thunderbolt shape) (**Figure 7.3**). The configuration of the ascending module delivery system is pre-curved (**Figure 7.2**). Both configurations allow accommodation to the native anatomy. The Descending Extension delivery system is straight.


Figure 7.1 – Delivery System Components

Figure 7.2 - Pre Curved Ascending Delivery System


Figure 7.3 - Pre Shaped Main Module Delivery System

3 | INDICATION FOR USE

The Nexus™ Aortic Arch Stent Graft System is indicated for the endovascular treatment of thoracic aortic diseases involving the aortic arch with proximal landing zone in ascending aorta and the Brachiocephalic artery. This includes:

- Aneurysm
- Dissecting aneurysm / dissection and intramural hematoma (IMH)
- False / Pseudo aneurysm if not infected
- Residual aneurysm/dissection following ascending aorta open repair
- Penetrating ulcer, if not infected

The Nexus™ Aortic Arch Stent Graft System is intended to exclude the lesion from the blood circulation in patients diagnosed with aortic arch/thoracic aortic diseases and who have appropriate anatomy to accommodate the Nexus™ system in an endovascular procedure.

4 | CONTRAINDICATIONS

- An acutely ruptured or unstable aneurysm.
- An acute vascular injury of the aorta due to trauma.
- An aneurysm that is mycotic or inflammatory.
- Connective tissue disease (e.g., Marfan's or Ehler's-Danlos syndromes).
- A "shaggy" aorta.
- Ascending aorta or Brachiocephalic artery that do not allow sufficient landing zone for the Nexus™ device.
- Patients with severe atherosclerosis or intraluminal thrombus of the aorta or the Brachiocephalic trunk.
- Patient whose arterial access sites (femoral and/or brachial/axillary) are not anticipated to accommodate the diameter of the Nexus™ Delivery System, due to size, tortuosity or hostile groins/arm(pit), unless conduit are used.
- Patient with a contraindication to undergo angiography.
- Patient with known sensitivities or allergies to the device materials.
- Patient who has a condition that threatens to infect the graft.

5 | WARNINGS AND PRECAUTIONS

Read all the instructions carefully. Failure to properly follow the instructions, warnings, and precautions may contribute to poor Nexus™ system performance which may lead to serious consequences and/or injuries to the patient, deterioration in patient's health and even patient's death.

5.1 General

- The Nexus™ Aortic Arch Stent Graft System should be used only by a physician's team with expertise in vascular interventional techniques, and trained in the use of this system. Specific training expectations are described in section 9.1, "Physician Training Requirements".
- Inadequate implantation of The Nexus™ Aortic Arch Stent Graft System may result in poor performance, severe complications and even death.
- A (cardio) vascular surgery team should always be available during endovascular grafting procedures in the event that conversion to open surgical repair is necessary.

5.2 Patient Selection

- Do not use the Nexus™ Aortic Arch Stent Graft System in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation procedure described in Section 10, DIRECTIONS FOR USE. Not following this instruction may lead to serious injuries to the patient, deterioration in patient's health and even patient's death.
- The Nexus™ Aortic Arch Stent Graft System is not recommended in patients who do not tolerate contrast agents necessary for intra-operative and post-operative follow-up imaging.
- **Key anatomic elements may affect successful exclusion of the thoracic lesion. These include:**
 - Ascending Aorta landing zone length (related to the chosen ascending module) ≥ 30 mm;
 - Ascending Aorta landing zone diameter should be: $29 \leq \varnothing \leq 39$ mm*
 - Descending Aorta landing zone length, either of the main module or the DE module: ≥ 30 mm
 - Descending Aorta landing zone diameter : $26 \leq \varnothing \leq 40$ mm

- Brachiocephalic Artery landing zone length (related to the chosen main module) ≥ 20 mm;
- Brachiocephalic Artery landing zone diameter: $11.5 \leq \varnothing \leq 18.5$ mm
- Take off angle between the Brachiocephalic Artery and the Aortic Arch perpendicular should be $\geq 125^\circ$, (see **Figure 8**).

*** For the shortest ascending module, ascending Aorta landing zone diameter should be: $29 \leq \varnothing \leq 36$ mm (please refer to Stent Graft Sizing Chart in section 2.2).**

- Additional key anatomic elements include: significant thrombus and/or calcium deposits at the implantation sites, specifically at the proximal aortic neck. Presence of calcification and/or plaque may compromise the fixation and sealing of the implantation sites. It may also promote emboli formation and subsequent adverse events. Necks that do not exhibit these key anatomic elements may be more conducive to graft migration.


Figure 8: The angle between the Brachiocephalic Artery and the Aortic arch

- In case of deployment of the Descending Extension (DE) Stent Graft, anatomical limitations that may affect successful treatment of the lesion include severe aortic angulation or if inter-componential overlap segment is completely unsupported by aortic wall. In the presence of anatomical limitations, a longer landing zone and additional stent grafts may be required to obtain adequate sealing and fixation.
- In case the difference between the proximal and distal intended landing diameters does not enable sufficient oversizing offered by a single DE endoprosthesis, a combination of multiple DE devices of different diameters is required.
- Inappropriate patient selection may lead to poor device performance or not in accordance with the specifications. This could result in injury to the patient, deterioration in patient's health and even patient's death.
- All patients should be advised that endovascular treatment requires lifelong, regular follow-up visits to assess patient's health and performance of the implanted endovascular stent graft. Patients with specific clinical findings (e.g. endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up.
- Patients experiencing reduced blood flow through the graft modules and/or leaks may be required to undergo secondary interventions or an open surgical procedure.

5.3 Before Implantation

- Pre-operative planning for vascular access, for Through & Through procedure and for device placement should be performed before the operation.
- A bolus of IV heparin should be administered before inserting the device as per standard of care to reduce the risk of thrombotic events.
- Carefully inspect the Nexus™ Aortic Arch Stent Graft System packaging and the delivery system for damage or defects prior to use. Do not use product if any sign of damage or breach of the sterile barrier is observed. Do not attempt to re-sterilize the Nexus™ Aortic Arch Stent Graft System or its components.
- Do not bend, kink, or otherwise challenge the Delivery System prior to implantation because it may cause deployment difficulties.

5.4 During Implantation

- Systemic anticoagulation should be used during the implantation procedure based on hospital or physician protocol.
- If heparin is contraindicated, an alternative anticoagulant should be considered.
- Exercise care in handling and in delivery technique to prevent blood vessel rupture.
- Blood vessels that are significantly calcified, stenosed, tortuous or thrombus-lined should be precluded from endovascular treatment.
- Do not advance the delivery system without having placed a guide wire first and without fluoroscopy. This may lead to major vascular complication, and even to patient's death.
- A suitable introducer sheath must be used for the Through-and-Through guide wire technique to accommodate the guide wire. In addition, the guide wire should remain in place and stressed until the Main Module is deployed in-situ in order to avoid migration of the main module during procedure. Migration can lead to inadequate sealing of the lesion and even to innominate artery occlusion which can eventually result in stroke and even death.
- Do not deploy the stent graft components in a location that can cause an endoleak or can occlude arteries necessary to supply blood flow to vital organs. This may necessitate surgical removal of the device.
- Use fluoroscopic guidance to advance the delivery system in order to detect kinking or alignment problems with the stent graft components. This may result in vascular injuries to the patient, deterioration in patient's health and even patient's death.
- Do not use excessive force to advance or withdraw the delivery system when resistance is encountered, it may lead to improper performance and may result in vascular injury to the patient and even patient's death. If resistance is felt during advancement of the guide wire or delivery system, stop and assess the cause of resistance. If the delivery system kinks during insertion, do not attempt to deploy the stent graft component. Remove the system and insert a new delivery system.
- Mild tension on the Through & Through wire can be applied during deployment of the Main Module to ensure it is straight and it does not fold; however, avoid excessive tension that might result in squeezing the gateway of the Ascending Module.
- Improper alignment of the radiopaque dot markers of the Ascending Module with the radiopaque large ring of the Main-Module and/or improper direction of "B" marker of the Ascending Module (must be turned toward the cranial side of the aorta) may impair the connection of the two modules and may result in endoleak or occlusion of the innominate artery which can result in stroke and even death.
- Improper alignment of the radiopaque "B" marker of the DE Module with the radiopaque dot marker of its overlapping module may impair the connection of the two modules and may result in endoleak or occlusion of the innominate artery which can result in stroke and even death.
- While in situ and for the purpose of correct orientation of the markers, do not rotate the Main Module and the Ascending Module delivery systems grip more than 360 degrees. In case the delivery system does not rotate, gently pull the delivery system backward out of the aortic arch and correct the marker orientation by rotating the delivery system grip (not more than 360 degrees), and advance the delivery system to the correct position. While in place, verify correct orientation before deployment.
- Inadequate seal zone may result in increased risk of leakage into the aneurysm in aneurysm treatment or migration of the stent graft.
- When deploying the stent graft, be sure to hold the front Grip of the delivery system stationary in order to ensure adequate implantation and avoid inaccurate positioning of the modules and as a result endoleak or occlusion of the innominate artery which can result in stroke and even death.
- Rapid Heart Pacing or other available techniques to reduce cardiac output and blood pressure, may be used during the deployment of the Ascending module.
- Ensure that the Main Module and the DE Module are deployed against the outer curve of the vasculature, e.g., by applying forward pressure on either the delivery system or both ends of its corresponding guidewire. As default, it is possible to deploy the overlapping segment of the DE module for primary fixation then carefully apply minor pressure on the DS to accommodate the outer curve before completing the deployment. **Note:** avoid excess force that might result in relative module migration or other hazardous disposition of adjacent implants.
- Stent graft components cannot be replaced or drawn back into the delivery system, even if the stent graft component is only partially deployed.
- If a balloon catheter is used, it should be inflated only inside the stent graft. Do not inflate outside the stent graft nor in the border area of the stent graft with the aortic wall. Follow the balloon catheter manufacturer Instructions for Use.
- Renal complications may occur from excess use of contrast agents, or as a result of embolization.
- Studies indicate that the danger of micro-embolization increases with increased procedure duration.^{1,2}
- Ensure that extra extensions and stents of suitable diameters are available in the operating room.

1) Becquemin JP, Majewski M, Fermani N , et al. Colon ischemia following abdominal aortic aneurysm repair in the era of endovascular abdominal aortic repair. J Vasc Surg 2008;47:258-63

2) Dadian N, Ohki T, Veith FJ, et al. Overt colon ischemia after endovascular aneurysm repair: the importance of microembolization as an etiology. J Vasc Surg. 2001;34(6):986-996.

5.5 Magnetic Resonance (MR) Imaging Safety

Non-clinical testing and analysis has demonstrated the Endospa[™] Nexus[™] Aortic Arch Stent-Graft System is MR conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 T and below
- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the Nexus[™] Aortic Arch Stent-Graft System is expected to produce a maximum temperature rise of less than 2 °C after 15 minutes of continuous scanning.

In non-clinical testing, the maximum image artifact caused by the device extends approximately 4 mm from the implant when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

The health state of the patient or the presence of other implants may require reduction of the MR limits.

6 | ANTICIPATED ADVERSE EVENTS

As for all medical devices and procedures, complications may occur with the implantation of the Nexus[™] Aortic Arch Stent Graft System. The following adverse events may occur for conventional similar procedures with endovascular stent-graft devices, as well as with the Nexus[™]. As the patient may undergo additional medical treatments for other underlying disease conditions, each of these medical conditions may cause a multitude of complications and adverse events. The potential adverse events, which may occur with the use of the Nexus[™] include, but are not limited to:

- 1) Acute respiratory distress syndrome (ARDS)
- 2) Allergies / inflammation - from residual EtO, contrast media-related toxicity
- 3) Amputation
- 4) Anesthetic complications and subsequent attendant problems
- 5) Aneurysm enlargement or rupture
- 6) Aortic damage, including perforation, dissection, bleeding, rupture, subsequent pericardial tamponade, aortic regurgitation, hemorrhagic pleural effusion, ischemia and death
- 7) Aortic valve damage and/or functional impairment
- 8) Vessel damage including perforation, dissection, bleeding, rupture subsequent ischemia, hemorrhage, thrombosis or pseudoaneurysm Obstruction/Occlusion of device or native vessel(s)
- 9) Arteriovenous fistula
- 10) Aortoesophageal fistula
- 11) Hemorrhage, requiring blood transfusion
- 12) Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension), Ventricular tachycardia or ventricular fibrillation
- 13) Claudication (e.g., buttock, lower limb)
- 14) Coronary arteries occlusion/stenosis
- 15) Death
- 16) Edema
- 17) Embolization (micro and macro) with transient or permanent ischemia or infarction
- 18) Emergency operation, Surgical conversion to open repair
- 19) Endoleak
- 20) Fever and localized / general inflammation
- 21) Peritonitis, Gastrointestinal complications (e.g., ileus, transient ischemia, infarction, necrosis)
- 22) Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection), Impotence
- 23) Infection of the aneurysm, device access site, including abscess formation, fever and pain, Sepsis
- 24) Ischemic effects to the limbs, kidneys, myocardium, mesentery, brain and spinal cord and subsequent tissue necrosis
- 25) Lymphatic complications and subsequent attendant problems (e.g., lymph fistula, lymphocele)
- 26) Multiple organ failure syndrome (MOF), Hepatic Failure
- 27) Myocardial perforation from a super-stiff guidewire
- 28) Neurologic complications, transient or permanent, local or systemic and subsequent attendant problems (e.g., Stroke, transient ischemic attack, syncope, lipothymy, paraplegia, paraparesis, paralysis, including vocal cord paralysis, diaphragmatic paralysis, Temporary Neurological Deficit (TND) and encephalopathy). Radiation exposure risks / repetitive radiation sessions
- 29) Pulmonary complications and subsequent attendant problems
- 30) Renal failure, renal complications and subsequent attendant problems (e.g., artery occlusion, contrast media toxicity, insufficiency)

- 31) Stent graft: improper component placement; incomplete component deployment; inadequate apposition, “beaking” and subsequent collapse of the component; clinically significant component migration; suture break; occlusion; infection; stent fracture; graft twisting or kinking; insertion and removal difficulties; graft material wear; dilatation; erosion; puncture and perigraft flow
- 32) Vascular access site complications, e.g. infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula, dissection, lymphocele, seroma.
- 33) Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- 34) Wound complications and subsequent attendant problems (e.g., dehiscence, infection, hematoma, seroma, cellulitis)

7 | TREATMENT AND FOLLOW-UP

7.1 Individualization of Treatment

Each Nexus™ Component must fit the patient’s anatomy and size. Proper sizing of the device is under the responsibility of the treating physician.

The aortic stent graft components should be over-sized approximately 10-20% from the inner diameter of the vessel or the Nexus™ stent graft component intended for extension. Please refer to section 2.2 for Device Size. The recommended overall length of multiple deployed in situ assembled modular components of the Nexus™ should extend from just below the bifurcation of Brachiocephalic-RCCA arteries, Aortic Arch to LCCA to Descending Aorta, and Aortic Arch to the Ascending Aorta, adjacently and distally to the Coronary Arteries. All lengths and diameters of the stent graft components necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are un-certain.

Physicians may consult Endospan Ltd. in their efforts to determine proper stent graft component dimensions based on the physician’s assessment of the patient’s anatomical measurements. The risks and benefits should be carefully considered for each patient before use of the Nexus™ Aortic Arch Stent Graft System.

Patient selection factors to be assessed should include:

- Patient’s age and life expectancy
- Families medical history
- Co-morbidities (e.g. cerebrovascular, cardiac, pulmonary or renal insufficiency prior to surgery)
- Patient morphologic suitability for endovascular repair

CAUTION: Vessel over-distension and damage may be caused by excessive oversizing of the stent graft in relation to the internal native blood vessel diameter.

CAUTION: Insufficient oversizing may increase risks of endoleaks and/or stent graft migration.

7.2 Treatment and Follow-up

- All patients with endovascular Aortic Arch or Thoracic Aorta repair require life-long, follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical finding (e.g., endoleaks, enlarging aneurysms or changes in structure or position of the endovascular graft) should receive enhanced follow-up.
- After endovascular graft placement, patients should have regularly follow-up visits which include periodic imaging to evaluate the stent graft, aneurysm size, and occlusion of vessels in the treatment area. In the case the Nexus™ is used to treat thoracic aneurysm, significant aneurysm enlargement (> 5 mm), appearance of a new endoleak, evidence of perigraft flow, or stent graft migration resulting in an inadequate seal zone should prompt initiate further investigation.
- Any endoleak left untreated during the implantation procedure must be carefully monitored after implantation.
- Additional treatment including endovascular treatment or surgical conversion should be strongly considered in the following cases:
 - Aneurysm growth > 5 mm (with or without endoleak) since last follow-up visit.
 - Persistent endoleak (with or without aneurysm growth).
 - Stent graft migration resulting in an inadequate seal zone.

7.3 Patient Counseling Information

The physician should review the following risks and benefits when counseling the patient about this endovascular device and procedure:

- Age and life expectancy.
- Risks and benefits related to open surgical repair.
- Risks and benefits related to endovascular repair.
- Differences between endovascular repair and surgical repair.
- Risks of spontaneous aneurysm rupture compared to endovascular repair.
- Risks related to non-interventional treatment (medical management).

- Possibility that subsequent endovascular or open surgical repair of the aneurysm may be required.
- Symptoms of aneurysm rupture.
- The long-term safety and effectiveness of aortic arch endovascular repair has not been established.
- Physicians should advise all patients that this treatment modality requires long-term, regular follow up to assess patients' health status and stent-graft performance.
- Physicians should provide description of symptoms of graft module occlusion, aneurysm enlargement or rupture
- Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should be monitored closely.
- Physicians must advise all patients that it is important to seek prompt medical attention if they experience signs of:
 - Symptoms of graft module occlusion, aneurysm enlargement or rupture.
 - Chest pain.
 - Chest pressure.
 - Difficulty in Breathing.
 - Tiredness.
 - Claudication.
 - Neurological symptoms as Confusion, syncope, stroke, transient ischemic attack, paraplegia, paraparesis, paralysis, including vocal cord paralysis.
- Patients should be counseled on the importance of adhering to the follow-up schedule.

8 | HOW SUPPLIED

8.1 Contents and Sterility

The Nexus™ Aortic Arch Stent Graft System includes two delivery systems (and additional delivery system in case of using the Descending Extension Module), each comprises an integral handle. Each delivery system is dedicated to and bears one of the stent graft modules. The Nexus™ is supplied sterilized by Ethylene Oxide (EtO), in a peel-open package.

The device is intended for single use only.

If the device is damaged or the integrity of the sterilization barrier has been compromised, do not use the product and contact Endospa Ltd. Do not reuse or attempt to re-sterilize.

8.2 Storage

Store the system at room temperature (between 15 °C to 25 °C), in a dark dry place.

9 | CLINICAL USE INFORMATION

9.1 Physician Training Requirements

All Physicians should complete training prior to using the Nexus™ Aortic Arch Stent Graft System.

CAUTION: The Nexus™ Aortic Arch Aneurysm Stent Graft System should be used only by physicians and teams trained in vascular and aortic interventional techniques and in the use of this device.

CAUTION: A (cardio) vascular surgery team should always be available during endovascular grafting procedures in the case a conversion to open surgical repair is necessary. It is strongly recommended that an interventional cardiologist is available in the case acute heart dysfunction as arrhythmia or pericardial hemorrhage occur.

Below are the skill/knowledge requirements for physicians using the Nexus™ System:

- Knowledge of the natural history of aortic arch and thoracic aortic diseases, in particular aortic arch aneurysms, and co-morbidities associated with aortic arch repair
- Radiographic, fluoroscopic, and angiographic image interpretation
- Appropriate use of radiographic contrast material
- Arterial cut-down, arteriotomy, and repair
- Percutaneous access and closure techniques
- Nonselective and selective guide wire and catheter techniques
- Embolization
- Angioplasty
- Endovascular (aortic) stent placement
- Techniques to minimize radiation exposure
- Device selection and sizing

- Rapid Heart Pacing, or other available techniques to reduce cardiac output and blood pressure

9.2 Inspection Prior Use

Inspect the device and packaging to verify that no damage or defect has occurred. If the “Use by” date has elapsed, if the device is damaged, or if the sterilization barrier has been compromised, do not use the device and contact Endospan Ltd. for return or replacement.

9.3 Additional Equipment Required

- Fluoroscopic imaging with digital angiographic capabilities (fixed units are recommended) and the ability to record and recall all imaging (including radiopaque ruler with centimeter) increments.
- Assorted guide wires and guiding catheters of adequate sizes.
- Assorted balloon catheters.
- Contrast media.
- Heparinized saline solution.
- Heart Pacing equipment, including electrodes (Cardiac Pacing)
- Large size stents of various diameters
- Stent graft extensions

9.4 Following Usage

After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy.

10 | DIRECTIONS FOR USE

The following instructions contain a basic guideline for device placement. These instructions are intended to guide the physician and do not take the place of physician judgment.

10.1 General Use Information

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and guide wires are applicable.

10.2 Pre-Implant Determinants

Correct sizing of the aortic arch, supra aortic and iliac arteries must be determined before implantation of the stent graft components using contrast-enhanced computer-aided tomography (CT) or MRA, as well as angiograms of vessels. 3D imaging may also be beneficial. These images should be available for review during the procedure. Vascular instruments and other surgical supplies needed to gain access to the artery should also be available.

Determine access paths and compatibility for an endovascular procedure through a 20 Fr sheath, the suitable device sizing based on patient’s anatomical measurements (see also Section 2.2 Device Size).

To reduce the risk of thromboembolism, it is recommended that the patient be heparinized for the duration of the procedure.

CAUTION: Never advance or retract equipment from the vasculature without using fluoroscopy. This may result in a major vascular complication and may lead to patient’s death.

10.3 Patient Preparation

- 10.3.1 Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
- 10.3.2 Position patient in supine position and imaging table allowing fluoroscopic visualization from the aortic arch to the iliac arteries.

10.4 Vascular Access

- 10.4.1 Establish bilateral femoral access (ipsilateral 7 Fr upgradable to 20 Fr, contralateral at least 5 Fr).
- 10.4.2 Establish right brachial/axillary access (7 Fr). Fixate the introducer sheath.

10.5 Device Preparation

Note: This section refers to both Main and Ascending Module delivery systems, as well as for the Descending Extension delivery system, if used.

- 10.5.1 Inspect the device packaging for any sign of damage or breach of the sterile barrier. If damage is observed, replace with another device.
- 10.5.2 Flush the delivery sheath with at least 50ml of heparinized saline through each of the three flushing ports (see **Figure 7.1** and **Figure 9 below**).

Note: For a longer and thorough flushing using a heparinized saline, connect the flushing device to the Stent graft flushing port, see figure below.

WARNING: Avoiding the flushing method proposed in this manual may result in an excess amount of air bubbles within the delivery system. This may result in major neurological complication, e.g., stroke.


Figure 9: Flushing through the stent graft flushing port

10.5.3 Prior to insertion, assess visualization and orientation of the radiopaque markers of the stent-graft under direct fluoroscopy of the delivery system. The radiopaque markers indicate the position of the proximal and distal edges of the graft material, and the orientation of modules.

10.5.4 Wet the shaft of the delivery system with heparinized saline to activate its hydrophilic coating prior or during the delivery system introduction. Materials Recommended for Through & Through procedure: Guide wire, stiff, 400/450 cm for through & through.

Note: Consumables are not included in the Nexus™ implantation system. Refer to Appendix A for the complete list of recommended material.

10.6 Through & Through

10.6.1 Perform a through & through; brachial-ipsilateral femoral, using guide wire 0.035", 400 cm (at least).

10.6.2 Fixate the wire in the brachial side, leaving ~30 cm emerging from the brachial introducer.

10.7 Main Module Insertion & Positioning

WARNING: Do not advance the delivery system without using a guide wire and fluoroscopy. This may result in a major vascular complication to the patient and even death.

Note: The vessels should be imaged as completely as possible, during the deployment of the Main Module, it is recommended to have two C-ARM configurations that allow accurate placement of the stent graft. (**Figure 10**)

Working view 1: "RCCA – RSA Plane" used for positioning of cranial end ("Brachiocephalic branch") of Main Module. This angle should allow the maximal angle between RCCA & RSA.

Working view 2: "BCT – ASC Plane" used for continued deployment of Main Module. This angle should allow the maximal angle between BCT & ASC.

Working View 1: RAO XX°, Caudal YY°

Definition - Maximize angle between RCC & RSA.

When in use - BCT segment deployment of Main Module.

**Working View 2: LAO XX°, Caudal YY°**

Definition - Maximize angle between BCT & ASC.

When in use - Completion deployment of Main Module.

Figure 10 – C-Arm Working Views (example)

- 10.7.1 Mount the delivery system over the through & through guide wire.
- 10.7.2 Just prior to insertion of the DS and during its insertion to the introducer sheath, perform simultaneous flushing with heparinized saline using a 20 ml syringe each, through the flushing ports according to the following stages: **(Figure 7.1)**
- Flush the DS flushing port until drops are seen in the stent graft port.
- Note:** Do not disconnect the injector at the end of this stage.
- Flush the SG flushing port simultaneously with flushing the DS flushing port, while advancing the DS into the introducer.
 - When reaching with the SG flushing port to introducer sheath remove SG flushing syringe.
- 10.7.3 Using continuous fluoroscopic guidance, advance the delivery system to the point the Fixation tube **(Figure 7.2)** emerges from the brachial/axillary introducer sheath. Ensure that the Main Module is deployed against the outer curve of the vasculature by applying forward pressure on either the delivery system or both ends of its corresponding guidewire.
- 10.7.4 Verify the radiopaque markers are visualized in the target location **(Figure 11)**:
- Caudal dot marker of the Brachiocephalic trunk is aligned with its ostia.
 - "B" marker is oriented with the cranial side of the Arch.
 - Verify the most cranial stent is below (caudal) the bifurcation of Brachiocephalic artery (bifurcation to RSA-RCC).
 - Radiopaque ring wire is facing the ascending arch.

Note: If the B marker is not positioned correctly rotate the delivery system grip no more than 360 degrees. In case the delivery system does not rotate, gently pull the delivery system backward out of the aortic arch and correct the marker orientation by rotating the delivery system grip (not more than 360 degrees). Then advance the delivery system to the correct position.

WARNING: Verify correct orientation of the markers before deployment. Failure to properly align the radiopaque markers may result in improper deployment of the stent graft.

Inject contrast media through an angiographic (pigtail) catheter into the Brachiocephalic/aortic arch and mark the position of the target location, either on the imaging screen.

- 10.7.5 Refine positioning of the stent graft so that the cranial edge of the graft fabric is just below the Brachiocephalic-RCCA bifurcation.

Caution: The angiographic catheter can be removed prior to deployment. However, if the angiographic catheter is not removed until after deployment, ensure that the catheter is not obstructing in any way the delivery or the stent graft.


Figure 11 – Radiopaque Markers of Main Module

10.8 Main Module Deployment

- 10.8.1 With one hand on the front grip, hold the delivery system stationary.
- 10.8.2 With the other hand, slowly withdraw the graft cover by rotating the deployment knob counter-clockwise, until the constrained cranial stent is exposed and 1 to 2 of the covered struts have been fully deployed (working view 1). **(Figure 12.a)**
- 10.8.3 Use angiography to verify position of the stent graft in relation to the Brachiocephalic-LCC bifurcation (working view 1).
- 10.8.4 Release cranial end of the Main-Module stent graft - capture mechanism by turning the capture release knob: One turn for releasing the capture at the cranial end of the Stent Graft.
- 10.8.5 Continue deployment of main module, allowing the self-projecting sleeve to open towards the ascending aorta, as proximally as possible/needed (working view 2) **(Figure 12.b)**.
- 10.8.6 Complete deployment by deploying the distal end of the stent graft either by rotating the deployment knob or by pressing the quick deployment button and pulling down (working view 2) **(Figure 13)**.

CAUTION: When using the quick deployment button to rapidly deploy the stent graft, be sure to hold the delivery system stationary in order to ensure adequate implantation and avoid inaccurate positioning of the modules and as a result endoleak or occlusion of the innominate artery which can result in stroke and even death.

CAUTION: Do not rotate the the delivery system after deployment had started as this may torque the device and cause it to rotate during deployment. This can cause a delivery system failure which will not allow to complete deployment and may result with procedure conversion.

CAUTION: If the graft cover is accidentally withdrawn, the stent graft will prematurely deploy and may be incorrectly positioned. This may lead to inaccurate positioning of the modules and as a result endoleak or occlusion of the innominate artery which can result in stroke and even death.

Figure 12.a | Cranial End Deployment (Working View 1)


Figure 12.b | Caudal Deployment (Working View 2)


Figure 12 – Deployment of Main Module


Figure 13 – Cranial Capture Release of Main Module

10.9 Main Module Delivery System Removal

CAUTION: In order to avoid module migration during this stage, verify, prior to delivery system removal and using fluoroscopy, that the Main module was completely deployed and there is no attachment between the module and the delivery system: **gently** try to attract the system, no movement of the module should be seen. Migration can lead to inadequate sealing of the lesion and in extreme cases to innominate artery occlusion which can eventually result in stroke and even death.

- 10.9.1 Gently apply tension and then apply fixation of the cranial end of the delivery system, "Fixation tube".
- 10.9.2 Withdraw the delivery system of the main module without releasing the fixation tube, until it is completely out of the main module.
- 10.9.3 Re-sheath the delivery system, distally to (below) the main module.
- 10.9.4 Remove completely the delivery system over the wire.
- 10.9.5 Release the through & through guide wire from the Iliac access and pull it partially back from the Brachial access, so that the soft extremity of the guide wire is located in the descending aorta.
- 10.9.6 Dilate balloon in the main module as needed.

10.10 Ascending Module Insertion & Positioning

- 10.10.1 Insert the delivery system over a guidewire and fluoroscopy.
- 10.10.2 Using continuous fluoroscopic guidance, advance the delivery system to Ascending aorta.
- 10.10.3 Adjust position of the C-arm so that the ring radiopaque wire is visualized as a line (**Figure 14**).
- 10.10.4 Verify the radiopaque markers are visualized in the target location (**Figure 14**):
 - a. Radiopaque dot markers are visualized in both ends of the graft.
 - b. "B" marker belly is oriented cranially.


Figure 14 – Radiopaque Markers of Ascending Module

- 10.10.5 Adjust position of the stent graft so that the distal end of the fabric, identified by the distal dot marker, is aligned with the radiopaque ring wire of the main module.
- 10.10.6 Inject contrast media through an angiographic (pigtail) catheter into the ascending aorta and mark the position of the target location, on the imaging screen.
- 10.10.7 Refine position of the stent graft so that the proximal end of the graft fabric is just distally to the Sino-Tubular Junction, verifying that the coronary arteries origin is not obstructed by the device. The distal end of the fabric, identified by a dot marker, is aligned with the radiopaque ring wire of the main module.

CAUTION: Please pay attention to floppy end of wire to reduce the risk of any damage or perforation to the heart and/or aorta, as traditionally addressed in similar procedures (e.g. TAVI). Alternatively, use techniques that provide stability and fixation of the guide wire inside the LV (such as snaring or looping).

10.11 Ascending Module Deployment

- 10.11.1 With one hand on the front grip, hold the delivery system stationary, while with the other hand slowly withdraw the graft cover by rotating the deployment knob counter-clockwise.
- 10.11.2 During unsheathing and module release, perform Rapid Heart Pacing.
 - Note:** Throughout the deployment verify alignment of the distal end of the fabric (identified by two dot markers) with the large tantalum oval marker of the main body (preferably in high magnification).
- 10.11.3 The proximal & distal capture mechanisms remain active after device deployment (**Figure 15**).


Figure 15a - Proximal capture release of Ascending Module


Figure 15b - Distal capture release of Ascending Module

Figure 15 – Capture Mechanism of Ascending Module

- 10.11.4 Release the capture mechanism by turning the capture release knob. Verify using fluoroscopic that both proximal and distal captures are released.

10.12 Ascending Module Delivery System Removal

CAUTION: In order to avoid module migration during this stage, verify, prior to delivery system removal and using fluoroscopy, that the module was completely deployed and there is no attachment between the module and the delivery system: **gently** try to retract the system, no movement of the module should be seen. Module migration can lead to inadequate sealing of the lesion and in extreme cases to innominate artery occlusion which can eventually result in stroke and even death.

- 10.12.1 Pull the delivery system until it is distally to the main module.
- 10.12.2 Re-sheath the delivery system and Remove the delivery system over the wire.
- 10.12.3 Procedure finalization: Dilate balloon as needed.

10.13 Descending Extension - Optional

- 10.13.1 Insert the delivery system over a 0.035" (0.89 mm) guidewire and using continuous fluoroscopic guidance, advance the delivery system to descending aorta.
- 10.13.2 Verify true entry into the already deployed module (either Main or DE) by applying slight tension on through-and-through wire or by using a super stiff wire.
- 10.13.3 For ideal visualization, position C-arm so that both distal dot markers of the already deployed module are aligned (relatively appear as one).
- 10.13.4 Acquire minimum overlap between the already deployed module and DE module by positioning the B marker of the DE in visual contact with the deployed module's distal dot marker.
Note: in Distal-to-Proximal order of deployment, the positioning is of the DE module's distal dot marker or Main Module's distal dot marker in alignment with the already deployed DE module's "B" marker (**Figure 16**).
- 10.13.5 With one hand on the front grip, hold the delivery system stationary, while with the other hand slowly withdraw the graft cover by rotating the deployment knob counter-clockwise.
- 10.13.6 Ensure that the DE module is deployed against the outer curve of the vasculature by applying forward pressure on either the delivery system or both ends of its corresponding guidewire. As default, it is possible to deploy the overlapping segment of the DE module for initial fixation then carefully apply minor pressure on the DS to accommodate the aortic outer curve and continue deployment.
Note: avoid excessive force that might result in module migration or other hazardous disposition of adjacent implants. Module migration can lead to endoleak and inadequate sealing of the lesion.
- 10.13.7 Once module's overlap segment is deployed and secured, move C-arm position to designated landing zone to follow and finalize the deployment with accuracy.
- 10.13.8 Release proximal end of the DE stent graft - capture mechanism by turning the capture release knob (arrow direction):
Two turns for releasing the proximal capture of the Stent Graft.

CAUTION: When using the quick deployment button to rapidly deploy the stent graft, be sure to hold the delivery system stationary in order to avoid inaccurate positioning of the module and as a result endoleak and inadequate sealing of the lesion.

CAUTION: Do not rotate the delivery system after deployment had started as this may torque the device and cause it to rotate during deployment. This can cause a delivery system failure which will not allow to complete deployment of the module and may result with procedure conversion.

CAUTION: If the graft cover is accidentally withdrawn, the stent graft will prematurely deploy and may be incorrectly positioned. This may lead to endoleak and inadequate sealing of the lesion.

- 10.13.9 Use fluoroscopic guidance to safely pull back the delivery system until the tip is in safe distance from the DE endoprosthesis.
- 10.13.10 Re-sheath the delivery system – this time, firmly hold the deployment knob stationary and simultaneously press the 'quick deployment button' while bringing the front grip backwards. Once completed, remove the delivery system out of the introducer and over the wire.


Figure 16 – Full assembly (Main Module, Ascending Module, and Descending Module)

10.14 Fixation Tube Removal

- 10.14.1 Pull the visible thread (Figure 17) until the knot is revealed, and cut the thread beyond the knot (cut only one suturing wire beyond the knot, not both).
- 10.14.2 Remove the fixation tube gently over the guide wire outside the introducer, until the thread is totally removed out of the introducer.


Figure 17 – Fixation Tube Release

11 | Bail-out Techniques

Stent Deployment: If five (5) counterclockwise rotations of the handle fail to initiate stent deployment (i.e. movement of the shaft tip and proximal stent strut deployment) conduct the following:

Press the release button and pull back the handle.

If after five centimeters, the shaft of the delivery system does not move backward (as seen under fluoroscopy): Remove the delivery system and replace it with a new one.

Release of the Ascending Module: In case of malfunction of the Ascending module release, make sure that the Through & Through wire is under tension and pull back the Ascending module into the main module until it is in the abdominal aorta

12 | FOLLOW-UP IMAGING RECOMMENDATION

Patients should be followed-up using conventional imaging techniques according to the institution's standard of care and physician discretion.

Current imaging of stent graft patients includes CT with a contrast medium. Alternative imaging modalities such as magnetic resonance imaging should be used in patients with impaired renal function or intolerance to contrast media.

Determination of imaging techniques should be based on the physician's clinical assessment of the patient and stent graft implant, including any adjunct procedures that may have been performed in conjunction with the endovascular stent graft procedure.

After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft.

Annual imaging is recommended to examine device integrity (stent fracture, separation between modules, to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease).

Contrast-enhanced CT should be used to assess stent graft fixation, deformation, apposition to the vessel wall at proximal and distal fixation sites, stent graft migration, stent graft patency, aneurysm size, occlusion of branch vessels, and endoleak (including source and type if present).

A pre-contrast scan of 1 mm thick slices is suggested to determine if there are calcifications or areas where metal artifacts may be misinterpreted as endoleak. An arterial phase with < 1 mm slice thickness and overlapping images is recommended. In aneurysms that are not shrinking and have no apparent endoleak or fixation problems, a late venous phase scan may be performed. The venous phase scan may also be performed with thicker collimation (5 mm). It is recommended that the source data set be archived in case specialized evaluation is needed later (volume measurements, 3-dimensional reconstruction, or computer-aided measurement software). If the aneurysm is not shrinking by more than 5 mm within the first year, volume measurements may be obtained as a more sensitive indicator of aneurysm size using 3-dimensional software.

The physician will determine the requirement pre-operative care for patients with allergies to contrast.

13 | Explanation of Symbols

| | | | |
|---|----------------------------------|---|---|
|  | Do Not Use If Package Is Damaged |  | Keep away from sunlight |
|  | Do Not Reuse |  | Keep Dry |
|  | CAUTION: see instruction for use |  | Storage Temperature Limitation |
| STERILE EO | Sterilized using Ethylene Oxide | EC REP | Authorized Representative in the European Community |
|  | Date of manufacture |  | Manufacturer |
|  | Use By |  | MR Conditional |
| SN | Serial Number |  | CE Marking of Conformity |
| REF | Catalogue Number | | |

14 | APPENDIX A – Consumables recommendation

The following products are recommended:

Guide Wires:

- 0.035 inch (0.89 mm) extra stiff wire guide, 300 cm; for example:
 - Cook Amplatz Ultra Stiff Wire Guides (AUS)
 - Cook Amplatz Extra Stiff Wire Guides (AES)
 - Cook Lunderquist Extra Stiff Wire Guides (LESDC)
- 0.035 inch Endo-Flex-teflon coated stainless steel 400 cm
- 0.035 inch (0.89 mm) standard wire guide; for example:
 - Cook Nimble™ Wire Guides
 - Cook Rosen curved Wire Guide
 - Terumo Radiofocus Guide wire-180, 260, 300, 400 cm
 - Terumo Angled Glidewire (hydrophilic, soft) 300 cm
- Extra support 0.014 coronary guide wire (e.g. Grand Slam) 300 cm

Molding Balloons:

- Molding Balloons (e.g., CODA) / 40 mm diameter
- POWERFLEX Balloons (e.g., Cordis)-12-16 mm

Introducer Sets – 18 Fr and 20 Fr; for example:

- Cook Check-Flo® Introducer Sets (18 Fr)
- Cook Extra Large Check-Flo Introducer Sets (20 Fr)
- Cook Flexor® Balkin Up & Over® Contralateral Introducers

Long Introducer; for example:

- Teleflex medical Europe I - arrow flex w/o g.w. 9 Fr x 100 cm
- Teleflex medical Europe I - arrow flex w/o g.w. 11 Fr x 80 cm

Sizing catheter; for example:

- Cook Aorous® Centimeter Sizing Catheters

Angiographic radiopaque tip catheters; for example:

- Cook Beacon® Tip Torcon NB Advantage Catheters (VERT) 125 cm
- Cook JR4 Torcon NB Advantage Catheters (JR4)
- Cook H1 Advantage Catheters (JR4)
- ST. JUDE Medical Fast cath 6^s Fr Flush Catheters
- Terumo angled glide catheter 5 Fr, 100 cm
- Marked Pig tail catheter 5 Fr
- Cross over catheter (e.g. omniflush 65 cm)
- Sheath 90 cm 9 Fr (e.g. Arrow)

Entry needles; for example:

- Cook single wall entry needles

Branch Stent – Examples:

- Atrium-Advanta V12 OTW (covered stent)

Bare metal stents - example:

- Palmaz

Additional aortic stent graft extensions

- Snares (Goose neck 30 mm, 10 mm)

RHP: Temporary transvenous pacemaker electrodes and a generator enable to generate up to 180 BPM.

15 | DISCLAIMER OF WARRANTY

ALTHOUGH THE NEXUS™ AORTIC ARCH STENT GRAFT SYSTEM HEREAFTER REFERRED TO AS THE 'PRODUCT', HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, ENDOSPAN LTD. AND THEIR RESPECTIVE AFFILIATES, (COLLECTIVELY "ENDOSPAN") HAVE NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. ENDOSPAN, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. ENDOSPAN SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND ENDOSPAN TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

THE EXCLUSIONS AND LIMITATIONS SET OUT ABOVE ARE NOT INTENDED TO, AND SHOULD NOT BE CONSTRUED SO AS TO CONTRAVENE MANDATORY PROVISIONS OF APPLICABLE LAW. IF ANY PART OR TERM OF THIS DISCLAIMER OF WARRANTY IS HELD TO BE ILLEGAL, UNENFORCEABLE OR IN CONFLICT WITH APPLICABLE LAW BY A COURT OF COMPETENT JURISDICTION, THE VALIDITY OF THE REMAINING PORTIONS OF THIS DISCLAIMER OF WARRANTY SHALL NOT BE AFFECTED, AND ALL RIGHTS AND OBLIGATIONS SHALL BE CONSTRUED AND ENFORCED AS IF THIS DISCLAIMER OF WARRANTY DID NOT CONTAIN THE PARTICULAR PART OR TERM HELD TO BE INVALID.

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MedNet GmbH,